

Data Sheet

Polygard[®]-CR Filters

Versatile depth filters for the clarification and prefiltration of process fluids

Polygard[®]-CR products are nominally rated filters designed for particle removal applications in liquids and gases. The graded-density depth structure of Polygard[®]-CR filter media provides maximum filtration capacity, and the all-polypropylene construction offers low extractables levels and broad chemical compatibility.

Benefits

- Broad range of micron ratings to match a wide array of particle removal applications
- Nominal particle retention rating
- Low extractables levels
- High capacity media
- Designed for rigorous process conditions and broad chemical compatibility
- Ideal for designing scalable solutions from bench top to full-scale manufacturing



Regulatory Compliance

Polygard®-CR filters are designed, developed, and manufactured in accordance with a Quality Management System approved by an accredited registering body to an ISO® 9000 Quality Systems Standard and are shipped with a Certificate of Quality. Each Opticap® XL capsule and cartridge filter is supported by a Validation Guide for compliance with regulatory requirements.

For traceability and easy identification, each filter is marked with identifying characteristics.

Multiple Formats Available

Polygard®-CR filters are available in two formats, eleven pore sizes, and multiple configurations that vary by filter sizes and the type of inlet and outlet connection.

Media Types

Polygard®-CR (nominal)

- 0.1 µm
- 0.3 µm
- 0.5 µm
- 1.0 µm
- 3.0 µm
- 5.0 µm
- 10.0 µm
- 25.0 µm
- 50.0 µm
- 75.0 µm
- 100.0 µm

Filter Formats

- Opticap® XL disposable capsule filters
- Cartridge filters



Opticap® XL Disposable Capsule Filters

Opticap® XL disposable capsule filters with Polygard®-CR media are available in multiple filtration sizes, providing an optimal choice for every application.

The patented Opticap® XL capsule design allows unparalleled thermal and hydraulic stress resistance in a disposable filter, resulting in reliability, high confidence in the sterility process and improved cleanliness. The unique capsule design with Polygard®-CR media minimizes hold-up volume and reduces production losses.

The Right Size

A wide range of filter sizes is available to fit all of your application needs for easy scale-up of your small volume filtration steps to larger, full-scale filtration processes.

The Right Connections

Self-contained and disposable, Opticap® XL capsule filters are supplied with a choice of inlet and outlet connections to optimize your filtration process, including sanitary flanges which enables a higher flow rate, fractional sanitary flanges, and hose barbs.



Cartridge Filters

Polygard®-CR cartridge filters provide high throughput and minimal differential pressure. Cartridges are robust, strong, resilient and are designed to withstand multiple steam-in-place cycles.

A full range of filter sizes is available to suit application requirements. A variety of connection options are offered for easy adaptation to existing housings.

Opticap® XL Filters




Cartridge Filters

Specifications


	Opticap® XL 1	Opticap® XL 5	Opticap® XL 10	Cartridge Filters 2-inch	Per 10-inch
Nominal dimensions					
Maximum length	21.6 cm (8.5 in.)	21.6 cm (8.5 in.)	33.5 cm (13.2 in.)	–	–
Diameter	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	10.7 cm (4.2 in.)	7.0 cm (2.8 in.)	7.0 cm (2.8 in.)
Filter element length	2.5 cm (1 in.)	12.7 cm (5 in.)	25.4 cm (10 in.)	–	–
Materials of construction filter media					
Filter media	Polypropylene			Polypropylene	
Supports	Polypropylene			Polypropylene	
Structural components	Polypropylene			Polypropylene	
Vent o-rings	Silicone			–	
O-rings	–			Silicone	
Cartridge code F gaskets	–			Polypropylene	
Vent/Drain	1/4 in. hose barb with double o-ring seal			–	
Maximum inlet pressure	5.5 bar (80 psi) at 25 °C			–	
	2.8 bar (40 psi) at 60 °C			–	
	1.0 bar (15 psi) at 80 °C			–	
Maximum operating temperature	–			80 °C continuous	
Maximum differential pressure					
Forward	4.8 bar (70 psid) at ambient room temperature			4.8 bar (70 psid) at 20 °C	
NVR gravimetric extractables	After autoclaving and a 24 hour soak in ASTM® Type 1 reagent grade water at controlled room temperature:				
	≤ 10 mg	≤ 30 mg	≤ 55 mg	–	≤ 50 mg
Bacterial endotoxin	Aqueous extraction contains < 0.5 EU/mL as determined by the Limulus Amebocyte Lysate (LAL) Test.				
Oxidizable substances	Capsules meet the requirements of the USP Oxidizable Substances Test after a water flush of:				
	≤ 1000 mL	≤ 2500 mL	≤ 5000 mL	–	5000 mL
Sterilization by autoclave	May be autoclaved for 3 cycles of 30 minutes at 126 °C. (Cannot be steam sterilized in-line.)			May be autoclaved for 10 cycles of 30 minutes at 126 °C or steam sterilized for 10 cycles for 30 minutes at 126 °C or hot water sanitized at 80 °C maximum for 30 minutes.	
Component material toxicity	Component materials were tested and meet the criteria of the USP <88> Reactivity Test for Class VI Plastics. Polygard®-CR filters meet the requirements of the USP <88> Safety Test utilizing a 0.9% sodium chloride extraction.				
Indirect food additive	All component materials meet the FDA Indirect Food Additive requirements cited in 21 CFR 177–182.				
European pressure Equipment directive	We certify that this product complies with the European Pressure Equipment Directive, 97/23/EC of 29 May 1997. This product has been classified under Article 3 § 3 of the Pressure Vessel Directive. It has been designed and manufactured in accordance with sound engineering practice to ensure safe use. In compliance with Article 3 § 3 of this Pressure Equipment Directive, this product does not bear the CE mark.			–	–

Typical Clean Water Flow Rates


Opticap® XL 1 Capsule with Polygard®-CR Media – 0.1 µm Nominal (KRK1)




Opticap® XL 5 Capsule with Polygard®-CR Media – 0.1 µm Nominal (KRK1)




Opticap® XL 10 Capsule with Polygard®-CR Media – 0.1 µm Nominal (KRK1)




Opticap® XL 1 Capsule with Polygard®-CR Media – 5.0 µm Nominal (KR05)



Opticap® XL 5 Capsule with Polygard®-CR Media – 5.0 µm Nominal (KR05)




Opticap® XL 10 Capsule with Polygard®-CR Media – 5.0 µm Nominal (KR05)



Opticap® XL Legends Refer to Connection Type

TT = 38 mm (1½ in.) Sanitary Flange Inlet and Outlet
 FF = 19 mm (¾ in.) Sanitary Flange Inlet and Outlet
 HH = 14 mm (9/16 in.) Hose Barb Inlet and Outlet

10-inch Cartridge Filters with Polygard®-CR Media




Cartridge Legend

Refers to Pore Size

CRK1 = 0.1 µm
 CRK3 = 0.3 µm
 CRA5 = 0.5 µm
 CR01 = 1.0 µm
 CR03 = 3.0 µm
 CR05 = 5.0 µm
 CR10 = 10.0 µm

Ordering Information

Opticap® XL Capsule Filters



*1-inch filter element in a 5-inch capsule housing.

Ordering Information

Cartridge Filters

<div style="border: 1px solid black; width: 30px; height: 30px; margin: 0 auto; display: flex; align-items: center; justify-content: center; font-weight: bold; font-size: 18px;">C</div> <p>Cartridge Filter</p>	<div style="border: 1px solid black; width: 30px; height: 30px; margin: 0 auto; display: flex; align-items: center; justify-content: center; font-weight: bold; font-size: 18px;">R</div> <p>Media Type R = Polygard®-CR</p>	<div style="display: flex; justify-content: space-around; width: 100%;"> <div style="border: 1px solid black; width: 30px; height: 30px; margin: 0 auto;"></div> <div style="border: 1px solid black; width: 30px; height: 30px; margin: 0 auto;"></div> </div> <p>Pore Size (nominal)</p> <p>K1 = 0.1 µm K3 = 0.3 µm A5 = 0.5 µm O1 = 1.0 µm O3 = 3.0 µm O5 = 5.0 µm O10 = 10.0 µm O25 = 25.0 µm O50 = 50.0 µm O75 = 75.0 µm O99 = 100.0 µm</p>	<div style="border: 1px solid black; width: 30px; height: 30px; margin: 0 auto;"></div> <p>Cartridge Code</p> <p>O = (2-222) O-ring 5 = (2-222) O-ring w/spear 7 = (2-226) O-ring w/locking tab and spear F = Double open end flat gasket M = (2-118) O-ring (2 in. cartridge only)</p>	<div style="border: 1px solid black; width: 30px; height: 30px; margin: 0 auto;"></div> <p>Cartridge Length</p> <p>1 = 10 inch 2 = 20 inch 3 = 30 inch 4 = 40 inch 0 = 2 inch (Cartridge Code M only)</p>	<div style="border: 1px solid black; width: 30px; height: 30px; margin: 0 auto;"></div> <p>O-ring Material</p> <p>O = Silicone (std) E = EP gasket V = Fluoroelastomers T = PTFE encapsulated fluoroelastomers P = Extended length version with polypropylene end caps and cage, silicone gaskets (Cartridge Code F only) C = Short length version with polypropylene end caps and cage, silicone gaskets (code F only) 2 = Silicone (Cartridge Code M only)</p>	<div style="display: flex; justify-content: space-around; width: 100%;"> <div style="border: 1px solid black; width: 30px; height: 30px; margin: 0 auto;"></div> <div style="border: 1px solid black; width: 30px; height: 30px; margin: 0 auto;"></div> </div> <p>Quantity per Package</p> <p>O3 = 3 per package (2 in. only) O6 = 6 per package (10, 20, 30 and 40 in.)</p>
---	--	--	---	---	--	---

*Not all configurations are available.



Merck Millipore and the M mark are trademarks of Merck KGaA, Darmstadt, Germany.
Polygard and Opticap are registered trademarks of Merck KGaA, Darmstadt, Germany.
ASTM is a trademark of American Society for Testing and Materials.
ISO is a registered trademark of the International Organization for Standardization.
Lit No. DS1838EN00 Rev. B 03/2011 Job No. DP SBU-11-05462 Printed in the U.S.A.
©2012 EMD Millipore Corporation, Billerica, MA 01821 USA. All rights reserved.

To Place an Order or Receive
Technical Assistance

lenntech

Rotterdamseweg 402m

2629HH Delft

The Netherlands

info@lenntech.com

www.lenntech.com

Tel. +31-15-261.09.00

Fax. +31-15-261.62.89

Polygard-CR-Filters

KRK1A51TT1	CRK101006	CRK301006	CRA501006	CR0101006	CR0301006	CR0501006	CR1001006	CR2501006	CR5001006	CR7501006	CR9901006	CRK1M0203
KRK1A51FF1	CRK101E06	CRK301E06	CRA501E06	CR0101E06	CR0301E06	CR0501E06	CR1001E06	CR2501E06	CR5001E06	CR7501E06	CR9901E06	CRK3M0203
KRK1A51HH1	CRK101V06	CRK301V06	CRA501V06	CR0101V06	CR0301V06	CR0501V06	CR1001V06	CR2501V06	CR5001V06	CR7501V06	CR9901V06	CRA5M0203
KRK1A05TT1	CRK101T06	CRK301T06	CRA501T06	CR0101T06	CR0301T06	CR0501T06	CR1001T06	CR2501T06	CR5001T06	CR7501T06	CR9901T06	CR01M0203
KRK1A05FF1	CRK102006	CRK302006	CRA502006	CR0102006	CR0302006	CR0502006	CR1002006	CR2502006	CR5002006	CR7502006	CR9902006	CR03M0203
KRK1A05HH1	CRK102E06	CRK302E06	CRA502E06	CR0102E06	CR0302E06	CR0502E06	CR1002E06	CR2502E06	CR5002E06	CR7502E06	CR9902E06	CR05M0203
KRK1A10TT1	CRK102V06	CRK302V06	CRA502V06	CR0102V06	CR0302V06	CR0502V06	CR1002V06	CR2502V06	CR5002V06	CR7502V06	CR9902V06	CR10M0203
KRK1A10FF1	CRK102T06	CRK302T06	CRA502T06	CR0102T06	CR0302T06	CR0502T06	CR1002T06	CR2502T06	CR5002T06	CR7502T06	CR9902T06	CR25M0203
KRK1A10HH1	CRK103006	CRK303006	CRA503006	CR0103006	CR0303006	CR0503006	CR1003006	CR2503006	CR5003006	CR7503006	CR9903006	CR50M0203
KRK3A51TT1	CRK103E06	CRK303E06	CRA503E06	CR0103E06	CR0303E06	CR0503E06	CR1003E06	CR2503E06	CR5003E06	CR7503E06	CR9903E06	CR75M0203
KRK3A51FF1	CRK103V06	CRK303V06	CRA503V06	CR0103V06	CR0303V06	CR0503V06	CR1003V06	CR2503V06	CR5003V06	CR7503V06	CR9903V06	CR99M0203
KRK3A51HH1	CRK103T06	CRK303T06	CRA503T06	CR0103T06	CR0303T06	CR0503T06	CR1003T06	CR2503T06	CR5003T06	CR7503T06	CR9903T06	
KRK3A05TT1	CRK104006	CRK304006	CRA504006	CR0104006	CR0304006	CR0504006	CR1004006	CR2504006	CR5004006	CR7504006	CR9904006	
KRK3A05FF1	CRK104E06	CRK304E06	CRA504E06	CR0104E06	CR0304E06	CR0504E06	CR1004E06	CR2504E06	CR5004E06	CR7504E06	CR9904E06	
KRK3A05HH1	CRK104V06	CRK304V06	CRA504V06	CR0104V06	CR0304V06	CR0504V06	CR1004V06	CR2504V06	CR5004V06	CR7504V06	CR9904V06	
KRK3A10TT1	CRK104T06	CRK304T06	CRA504T06	CR0104T06	CR0304T06	CR0504T06	CR1004T06	CR2504T06	CR5004T06	CR7504T06	CR9904T06	
KRK3A10FF1	CRK151006	CRK351006	CRA551006	CR0151006	CR0351006	CR0551006	CR1051006	CR2551006	CR5051006	CR7551006	CR9951006	
KRK3A10HH1	CRK151E06	CRK351E06	CRA551E06	CR0151E06	CR0351E06	CR0551E06	CR1051E06	CR2551E06	CR5051E06	CR7551E06	CR9951E06	
KRA5A51TT1	CRK151V06	CRK351V06	CRA551V06	CR0151V06	CR0351V06	CR0551V06	CR1051V06	CR2551V06	CR5051V06	CR7551V06	CR9951V06	
KRA5A51FF1	CRK151T06	CRK351T06	CRA551T06	CR0151T06	CR0351T06	CR0551T06	CR1051T06	CR2551T06	CR5051T06	CR7551T06	CR9951T06	
KRA5A51HH1	CRK152006	CRK352006	CRA552006	CR0152006	CR0352006	CR0552006	CR1052006	CR2552006	CR5052006	CR7552006	CR9952006	
KRA5A05TT1	CRK152E06	CRK352E06	CRA552E06	CR0152E06	CR0352E06	CR0552E06	CR1052E06	CR2552E06	CR5052E06	CR7552E06	CR9952E06	
KRA5A05FF1	CRK152V06	CRK352V06	CRA552V06	CR0152V06	CR0352V06	CR0552V06	CR1052V06	CR2552V06	CR5052V06	CR7552V06	CR9952V06	
KRA5A05HH1	CRK152T06	CRK352T06	CRA552T06	CR0152T06	CR0352T06	CR0552T06	CR1052T06	CR2552T06	CR5052T06	CR7552T06	CR9952T06	
KRA5A10TT1	CRK153006	CRK353006	CRA553006	CR0153006	CR0353006	CR0553006	CR1053006	CR2553006	CR5053006	CR7553006	CR9953006	
KRA5A10FF1	CRK153E06	CRK353E06	CRA553E06	CR0153E06	CR0353E06	CR0553E06	CR1053E06	CR2553E06	CR5053E06	CR7553E06	CR9953E06	
KRA5A10HH1	CRK153V06	CRK353V06	CRA553V06	CR0153V06	CR0353V06	CR0553V06	CR1053V06	CR2553V06	CR5053V06	CR7553V06	CR9953V06	
KR01A51TT1	CRK153T06	CRK353T06	CRA553T06	CR0153T06	CR0353T06	CR0553T06	CR1053T06	CR2553T06	CR5053T06	CR7553T06	CR9953T06	
KR01A51FF1	CRK154006	CRK354006	CRA554006	CR0154006	CR0354006	CR0554006	CR1054006	CR2554006	CR5054006	CR7554006	CR9954006	
KR01A51HH1	CRK154E06	CRK354E06	CRA554E06	CR0154E06	CR0354E06	CR0554E06	CR1054E06	CR2554E06	CR5054E06	CR7554E06	CR9954E06	
KR01A05TT1	CRK154V06	CRK354V06	CRA554V06	CR0154V06	CR0354V06	CR0554V06	CR1054V06	CR2554V06	CR5054V06	CR7554V06	CR9954V06	
KR01A05FF1	CRK154T06	CRK354T06	CRA554T06	CR0154T06	CR0354T06	CR0554T06	CR1054T06	CR2554T06	CR5054T06	CR7554T06	CR9954T06	
KR01A05HH1	CRK171006	CRK371006	CRA571006	CR0171006	CR0371006	CR0571006	CR1071006	CR2571006	CR5071006	CR7571006	CR9971006	
KR01A10TT1	CRK171E06	CRK371E06	CRA571E06	CR0171E06	CR0371E06	CR0571E06	CR1071E06	CR2571E06	CR5071E06	CR7571E06	CR9971E06	
KR01A10FF1	CRK171V06	CRK371V06	CRA571V06	CR0171V06	CR0371V06	CR0571V06	CR1071V06	CR2571V06	CR5071V06	CR7571V06	CR9971V06	
KR01A10HH1	CRK171T06	CRK371T06	CRA571T06	CR0171T06	CR0371T06	CR0571T06	CR1071T06	CR2571T06	CR5071T06	CR7571T06	CR9971T06	
KR03A51TT1	CRK172006	CRK372006	CRA572006	CR0172006	CR0372006	CR0572006	CR1072006	CR2572006	CR5072006	CR7572006	CR9972006	
KR03A51FF1	CRK172E06	CRK372E06	CRA572E06	CR0172E06	CR0372E06	CR0572E06	CR1072E06	CR2572E06	CR5072E06	CR7572E06	CR9972E06	
KR03A51HH1	CRK172V06	CRK372V06	CRA572V06	CR0172V06	CR0372V06	CR0572V06	CR1072V06	CR2572V06	CR5072V06	CR7572V06	CR9972V06	
KR03A05TT1	CRK172T06	CRK372T06	CRA572T06	CR0172T06	CR0372T06	CR0572T06	CR1072T06	CR2572T06	CR5072T06	CR7572T06	CR9972T06	
KR03A05FF1	CRK173006	CRK373006	CRA573006	CR0173006	CR0373006	CR0573006	CR1073006	CR2573006	CR5073006	CR7573006	CR9973006	
KR03A05HH1	CRK173E06	CRK373E06	CRA573E06	CR0173E06	CR0373E06	CR0573E06	CR1073E06	CR2573E06	CR5073E06	CR7573E06	CR9973E06	
KR03A10TT1	CRK173V06	CRK373V06	CRA573V06	CR0173V06	CR0373V06	CR0573V06	CR1073V06	CR2573V06	CR5073V06	CR7573V06	CR9973V06	
KR03A10FF1	CRK173T06	CRK373T06	CRA573T06	CR0173T06	CR0373T06	CR0573T06	CR1073T06	CR2573T06	CR5073T06	CR7573T06	CR9973T06	
KR03A10HH1	CRK174006	CRK374006	CRA574006	CR0174006	CR0374006	CR0574006	CR1074006	CR2574006	CR5074006	CR7574006	CR9974006	
KR05A51TT1	CRK174E06	CRK374E06	CRA574E06	CR0174E06	CR0374E06	CR0574E06	CR1074E06	CR2574E06	CR5074E06	CR7574E06	CR9974E06	
KR05A51FF1	CRK174V06	CRK374V06	CRA574V06	CR0174V06	CR0374V06	CR0574V06	CR1074V06	CR2574V06	CR5074V06	CR7574V06	CR9974V06	
KR05A51HH1	CRK174T06	CRK374T06	CRA574T06	CR0174T06	CR0374T06	CR0574T06	CR1074T06	CR2574T06	CR5074T06	CR7574T06	CR9974T06	
KR05A05TT1	CRK1F1P06	CRK3F1P06	CRA5F1P06	CR01F1P06	CR03F1P06	CR05F1P06	CR10F1P06	CR25F1P06	CR50F1P06	CR75F1P06	CR99F1P06	
KR05A05FF1	CRK1F2P06	CRK3F2P06	CRA5F2P06	CR01F2P06	CR03F2P06	CR05F2P06	CR10F2P06	CR25F2P06	CR50F2P06	CR75F2P06	CR99F2P06	
KR05A05HH1	CRK1F3P06	CRK3F3P06	CRA5F3P06	CR01F3P06	CR03F3P06	CR05F3P06	CR10F3P06	CR25F3P06	CR50F3P06	CR75F3P06	CR99F3P06	
KR05A10TT1	CRK1F4P06	CRK3F4P06	CRA5F4P06	CR01F4P06	CR03F4P06	CR05F4P06	CR10F4P06	CR25F4P06	CR50F4P06	CR75F4P06	CR99F4P06	
KR05A10FF1	CRK1F1C06	CRK3F1C06	CRA5F1C06	CR01F1C06	CR03F1C06	CR05F1C06	CR10F1C06	CR25F1C06	CR50F1C06	CR75F1C06	CR99F1C06	
KR05A10HH1	CRK1F2C06	CRK3F2C06	CRA5F2C06	CR01F2C06	CR03F2C06	CR05F2C06	CR10F2C06	CR25F2C06	CR50F2C06	CR75F2C06	CR99F2C06	
KR10A51TT1	CRK1F3C06	CRK3F3C06	CRA5F3C06	CR01F3C06	CR03F3C06	CR05F3C06	CR10F3C06	CR25F3C06	CR50F3C06	CR75F3C06	CR99F3C06	
KR10A51FF1	CRK1F4C06	CRK3F4C06	CRA5F4C06	CR01F4C06	CR03F4C06	CR05F4C06	CR10F4C06	CR25F4C06	CR50F4C06	CR75F4C06	CR99F4C06	
KR10A51HH1	KR75A10HH1	KR75A51FF1	KR50A05HH1	KR50A51TT1	KR25A05FF1	KR99A10TT1	KR25A51FF1	KR10A10FF1	KR75A10TT1	KR99A05FF1	KR25A10FF1	
KR10A05TT1	KR99A51TT1	KR75A51HH1	KR50A10TT1	KR50A51FF1	KR25A05HH1	KR99A10FF1	KR25A51HH1	KR10A10HH1	KR75A10FF1	KR99A05HH1	KR25A10HH1	
KR10A05FF1	KR99A51FF1	KR75A05TT1	KR50A10FF1	KR50A51HH1	KR25A10TT1	KR99A10HH1	KR25A05TT1	KR25A51TT1	KR75A05HH1	KR99A05TT1	KR50A05FF1	
KR10A05HH1	KR99A51HH1	KR75A05FF1	KR50A10HH1	KR50A05TT1	KR10A10TT1	KR99A05TT1						

Lenntech B.V.
T +31-15-261.09.00
F +31-15-261.62.89
info@lenntech.com
www.lenntech.com